

# Bruker NMR Cheat Sheet

## To start automation:

- 1) Open the Topspin program
- 2) Type "iconnmr" into the text entry box that runs across the bottom of the Topspin window.
- 3) Select "Automation" then select your user account from the pop-up list
- 4) Double click on the line that corresponds to your sample position in the autosampler.
- 5) After filling out all of the sample information, make sure your line is highlighted then click "submit"
- 6) At the top of the window, click start
- 7) Enter the number that corresponds to your sample position in the autosampler, then click "ok".

## Commands to run a sample manually from TopSpin:

- 1) "sx #" to inject your sample into the magnet ("#" should be the number that corresponds to your sample's position in the autosampler, so, for example; "sx 3" would inject sample 3).
- 2) "edc" to enter the requisite sample information, select the experiment, and create a data set.
- 3) "rsh" then select "LASTBEST" to read the most recent good shim file into memory.
- 4) "lock" then select your solvent to lock.
- 5) "atma" on the 400 or 600 to run the automatic tune program (it is a fixed tune on the 300).
- 6) "getprosol" to read the probehead and solvent-dependent parameters
- 7) "topshim" to run automated shimming
- 8) "rga" to automatically set the receiver gain
- 9) "zg" to start an acquisition
- 10) "ft" then "apk" for a quick look at your spectrum
- 11) "sx ej" to retrieve your sample

## *Other Commands that could be useful:*

- "eda" to edit data acquisition parameters
- "tr" to transfer data to disk during an acquisition
- "rpar" to read a parameter set
- "abs" automatic baseline correction
- "em" exponential window multiplication of the FID
- "efp" exponential window multiplication followed by a Fourier transform and phase correction
- "gm" Gaussian window multiplication (gmp for gm + ft + apk)
- "ppp" to automatically pick peaks